



Caprin1 Antibody

Product Code	CSB-PA711733HA01RA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q5M9G3
Immunogen	Recombinant Rat Caprin-1 protein (2-707AA)
Raised In	Rabbit
Species Reactivity	Rat
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Caprin-1 (Cytoplasmic activation- and proliferation-associated protein 1) (GPI-anchored protein p137) (GPI-p137) (p137GPI) (RNA granule protein 105), Caprin1, Gpiap1 Rng105
Immunogen Species	Rattus norvegicus (Rat)
Research Area	Cell Biology
Target Names	Caprin1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.