





# FOXP2 Antibody

<b>Product Code</b>	CSB-PA008842DSR1HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	O15409
<b>Immunogen</b>	Recombinant Human Forkhead box protein P2 protein (416-715AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Tested Applications</b>	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
<b>Purification Method</b>	Antigen Affinity Purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Forkhead box protein P2 (CAG repeat protein 44) (Trinucleotide repeat-containing gene 10 protein), FOXP2, CAGH44 TNRC10
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Research Area</b>	Epigenetics and Nuclear Signaling
<b>Target Names</b>	FOXP2

## Image



Immunohistochemistry of paraffin-embedded human prostate tissue using CSB-PA008842DSR1HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human placenta tissue using CSB-PA008842DSR1HU at dilution of 1:100