





FOXP2 Antibody

Product Code	CSB-PA008842DSR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O15409
Immunogen	Recombinant Human Forkhead box protein P2 protein (416-715AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Forkhead box protein P2 (CAG repeat protein 44) (Trinucleotide repeat-containing gene 10 protein), FOXP2, CAGH44 TNRC10
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	FOXP2

Image



Immunohistochemistry of paraffin-embedded human prostate tissue using CSB-PA008842DSR1HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human placenta tissue using CSB-PA008842DSR1HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.