



# acpP Antibody, HRP conjugated

Product Code	CSB-PA015636DB01SKX
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P0A001
Immunogen	Recombinant Staphylococcus aureus Acyl carrier protein (1-77AA)
Raised In	Rabbit
Species Reactivity	Staphylococcus aureus
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Acyl carrier protein (ACP), acpP, hmrB
Immunogen Species	Staphylococcus aureus (strain Mu50 / ATCC 700699)
Research Area	Others
Target Names	acpP