



Gad2 Antibody

Product Code	CSB-PA009160HA01MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P48320
Immunogen	Recombinant Mouse Glutamate decarboxylase 2 protein (1-585AA)
Raised In	Rabbit
Species Reactivity	Mouse
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Glutamate decarboxylase 2 (EC 4.1.1.15) (65 kDa glutamic acid decarboxylase) (GAD-65) (Glutamate decarboxylase 65 kDa isoform), Gad2, Gad65
Immunogen Species	Mus musculus (Mouse)
Research Area	Neuroscience
Target Names	Gad2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.