



Gad2 Antibody

Product Code	CSB-PA009160HA01MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P48320
Immunogen	Recombinant Mouse Glutamate decarboxylase 2 protein (1-585AA)
Raised In	Rabbit
Species Reactivity	Mouse
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Glutamate decarboxylase 2 (EC 4.1.1.15) (65 kDa glutamic acid decarboxylase) (GAD-65) (Glutamate decarboxylase 65 kDa isoform), Gad2, Gad65
Immunogen Species	Mus musculus (Mouse)
Research Area	Neuroscience
Target Names	Gad2