



THAP11 Antibody

Product Code	CSB-PA846612LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96EK4
Immunogen	Recombinant Human THAP domain-containing protein 11 protein (1-313AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	THAP domain-containing protein 11, THAP11
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	THAP11
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.