




WDR77 Antibody

Product Code	CSB-PA859030LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BQA1
Immunogen	Recombinant Human Methylosome protein 50 protein (1-312AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Methylosome protein 50 (MEP-50) (Androgen receptor cofactor p44) (WD repeat-containing protein 77) (p44/Mep50), WDR77, MEP50 WD45
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	WDR77

Image



Immunohistochemistry of paraffin-embedded human skeletal muscle tissue using CSB-PA859030LA01HU at dilution of 1:100

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.