



TGFB1 Antibody

Product Code	CSB-PA023446LA01CH
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P09531
Immunogen	Recombinant Chicken Transforming growth factor beta-1 proprotein protein (260-373AA)
Raised In	Rabbit
Species Reactivity	Chicken
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Transforming growth factor beta-1 proprotein (TGF-beta-1) [Cleaved into: Latency-associated peptide (LAP)] (Transforming growth factor beta-1), TGFB1
Immunogen Species	Chicken
Research Area	Others
Target Names	TGFB1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.