



CELA3A Antibody, Biotin conjugated

Product Code	CSB-PA007588LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P09093
Immunogen	Recombinant Human Chymotrypsin-like elastase family member 3A protein (29-270AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Chymotrypsin-like elastase family member 3A (EC 3.4.21.70) (Elastase IIIA) (Elastase-3A) (Protease E), CELA3A, ELA3 ELA3A
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	CELA3A
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.