




SIRPD Antibody

Product Code	CSB-PA021337LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H106
Immunogen	Recombinant Human Signal-regulatory protein delta protein (30-197AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB; Recommended dilution: WB:1:1000-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Signal-regulatory protein delta (SIRP-delta) (Protein tyrosine phosphatase non-receptor type substrate 1-like 2), SIRPD, PTPNS1L2
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	SIRPD

Image



Western blot
 All lanes: SIRPD antibody at 14µg/ml
 Lane 1: 293T whole cell lysate
 Lane 2: Mouse spleen tissue
 Secondary
 Goat polyclonal to rabbit IgG at 1/10000 dilution
 Predicted band size: 22 kDa
 Observed band size: 22 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.