



nikR Antibody

Product Code	CSB-PA363929HA01EOD
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P0A6Z6
Immunogen	Recombinant Escherichia coli Nickel-responsive regulator protein (1-133AA)
Raised In	Rabbit
Species Reactivity	Escherichia coli
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Nickel-responsive regulator, nikR, yhhG
Immunogen Species	Escherichia coli (strain K12)
Research Area	Others
Target Names	nikR