



nikR Antibody, FITC conjugated

| | |
|----------------------------|--|
| Product Code | CSB-PA363929HC01EOD |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P0A6Z6 |
| Immunogen | Recombinant Escherichia coli Nickel-responsive regulator protein (1-133AA) |
| Raised In | Rabbit |
| Species Reactivity | Escherichia coli |
| Form | Liquid |
| Conjugate | FITC |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Nickel-responsive regulator, nikR, yhhG |
| Immunogen Species | Escherichia coli (strain K12) |
| Research Area | Others |
| Target Names | nikR |