



CD2 Antibody, HRP conjugated

Product Code	CSB-PA004894LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P06729
Immunogen	Recombinant Human T-cell surface antigen CD2 protein (25-209AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	T-cell surface antigen CD2 (Erythrocyte receptor) (LFA-2) (LFA-3 receptor) (Rosette receptor) (T-cell surface antigen T11/Leu-5) (CD antigen CD2), CD2, SRBC
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	CD2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.