



Gpi Antibody, Biotin conjugated

Product Code	CSB-PA009717LD01MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P06745
Immunogen	Recombinant Mouse Glucose-6-phosphate isomerase protein (2-558AA)
Raised In	Rabbit
Species Reactivity	Mouse
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Glucose-6-phosphate isomerase (GPI) (EC 5.3.1.9) (Autocrine motility factor) (AMF) (Neuroleukin) (NLK) (Phosphoglucose isomerase) (PGI) (Phosphohexose isomerase) (PHI), Gpi, Gpi1
Immunogen Species	Mus musculus (Mouse)
Research Area	Neuroscience
Target Names	Gpi
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.