



# TES-26 Antibody

|                     |  |
|---------------------|--|
| Product Code        | CSB-PA347478HA01THA  |
| Storage             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.                      |
| Uniprot No.         | P54190   |
| Immunogen           | Recombinant Toxocara canis 26 kDa secreted antigen protein (22-262AA)              |
| Raised In           | Rabbit   |
| Species Reactivity  | Toxocara canis   |
| Tested Applications | ELISA  |
| Form                | Liquid   |
| Conjugate           | Non-conjugated   |
| Storage Buffer      | Preservative: 0.03% Proclin 300<br>Constituents: 50% Glycerol, 0.01M PBS, PH 7.4   |
| Purification Method | >95%, Protein G purified   |
| Isotype             | IgG  |
| Clonality           | Polyclonal   |
| Alias               | 26 kDa secreted antigen (Toxocara excretory-secretory antigen 26) (TES-26), TES-26 |
| Immunogen Species   | Toxocara canis (Canine roundworm)  |
| Research Area       | Others   |
| Target Names        | TES-26   |