



TES-26 Antibody

Product Code	CSB-PA347478HA01THA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P54190
Immunogen	Recombinant Toxocara canis 26 kDa secreted antigen protein (22-262AA)
Raised In	Rabbit
Species Reactivity	Toxocara canis
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	26 kDa secreted antigen (Toxocara excretory-secretory antigen 26) (TES-26), TES-26
Immunogen Species	Toxocara canis (Canine roundworm)
Research Area	Others
Target Names	TES-26
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.