



# CAPN2 Antibody

<b>Product Code</b>	CSB-PA004496LA01PI
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P43367
<b>Immunogen</b>	Recombinant Pig Calpain-2 catalytic subunit protein (1-324AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Pig
<b>Tested Applications</b>	ELISA
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Calpain-2 catalytic subunit (EC 3.4.22.53) (Calcium-activated neutral proteinase 2) (CANP 2) (Calpain M-type) (Calpain-2 large subunit) (Millimolar-calpain) (M-calpain) (Fragment), CAPN2
<b>Immunogen Species</b>	Pig
<b>Research Area</b>	Others
<b>Target Names</b>	CAPN2
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.