



vpx Antibody

| | |
|----------------------------|-----------------------------------------------------------------------------------------|
| Product Code | CSB-PA325553LA01HLU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P18099 |
| Immunogen | Recombinant Human immunodeficiency virus type 2 subtype A Protein Vpx protein (1-113AA) |
| Raised In | Rabbit |
| Species Reactivity | Human immunodeficiency virus type 2 subtype A |
| Tested Applications | ELISA |
| Form | Liquid |
| Conjugate | Non-conjugated |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Protein Vpx (Viral protein X) (X ORF protein), vpx |
| Immunogen Species | Human immunodeficiency virus type 2 subtype A (isolate BEN) (HIV-2) |
| Research Area | Others |
| Target Names | vpx |