



vpx Antibody, HRP conjugated

Product Code	CSB-PA325553LB01HLU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P18099
Immunogen	Recombinant Human immunodeficiency virus type 2 subtype A Protein Vpx protein (1-113AA)
Raised In	Rabbit
Species Reactivity	Human immunodeficiency virus type 2 subtype A
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein Vpx (Viral protein X) (X ORF protein), vpx
Immunogen Species	Human immunodeficiency virus type 2 subtype A (isolate BEN) (HIV-2)
Research Area	Others
Target Names	vpx
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.