



SFTA3 Antibody

Product Code	CSB-PA021167LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P0C7M3
Immunogen	Recombinant Human Surfactant-associated protein 3 protein (1-94AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Surfactant-associated protein 3 (Surfactant-associated protein H) (SP-H), SFTA3, SFTPH
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	SFTA3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.