



atpE Antibody

Product Code	CSB-PA428118LA01AUQ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	A8MJW4
Immunogen	Recombinant strain OhILAs ATP synthase subunit c protein (1-88AA)
Raised In	Rabbit
Species Reactivity	strain OhILAs
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	ATP synthase subunit c (ATP synthase F(0) sector subunit c) (F-type ATPase subunit c) (F-ATPase subunit c) (Lipid-binding protein), atpE
Immunogen Species	Alkaliphilus oremlandii (strain OhILAs) (Clostridium oremlandii (strain OhILAs))
Research Area	Others
Target Names	atpE
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.