



MOB1A Antibody, FITC conjugated

Product Code	CSB-PA872519HC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H8S9
Immunogen	Recombinant Human MOB kinase activator 1A protein (2-216AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	MOB kinase activator 1A (Mob1 alpha) (Mob1A) (Mob1 homolog 1B) (Mps one binder kinase activator-like 1B), MOB1A, C2orf6 MOB4B MOBK1B MOBKL1B
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	MOB1A