



Omega-gliadin Antibody, HRP conjugated

Product Code	CSB-PA360898DB01TKR
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P02865
Immunogen	Recombinant Triticum monococcum Omega-gliadin protein (1-28AA)
Raised In	Rabbit
Species Reactivity	Triticum monococcum
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Omega-gliadin (Fragment)
Immunogen Species	Triticum monococcum (Einkorn wheat) (Crithodium monococcum)
Research Area	Others
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.