



Reg3g Antibody, Biotin conjugated

Product Code	CSB-PA019549LD01RA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P42854
Immunogen	Recombinant Rat Regenerating islet-derived protein 3-gamma protein (27-174AA)
Raised In	Rabbit
Species Reactivity	Rat
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Regenerating islet-derived protein 3-gamma (REG-3-gamma) (Pancreatitis-associated protein 3) (Regenerating islet-derived protein III-gamma) (Reg III-gamma) [Cleaved into: Regenerating islet-derived protein 3-gamma 16.5 kDa form; Regenerating islet-derived protein 3-gamma 15 kDa form], Reg3g, Pap3
Immunogen Species	Rattus norvegicus (Rat)
Research Area	Immunology
Target Names	Reg3g