



ORF2 Antibody, HRP conjugated

Product Code	CSB-PA727851LB01HVZ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q68985
Immunogen	Recombinant Hepatitis E virus genotype 1 Secreted protein ORF2 protein (20-660AA)
Raised In	Rabbit
Species Reactivity	Hepatitis E virus genotype 1
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Secreted protein ORF2 (Protein ORF2) (pORF2) ORF2
Immunogen Species	Hepatitis E virus genotype 1
Research Area	Others
Target Names	ORF2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.