




INHBA Antibody

Product Code	CSB-PA07589A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P08476
Immunogen	Recombinant Human Inhibin beta A chain protein (311-426AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:500-1:5000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Inhibin beta A chain (Activin beta-A chain) (Erythroid differentiation protein) (EDF), INHBA
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	INHBA

Image



Western blot

All lanes: INHBA antibody at 2µg/ml

Lane 1: Hela whole cell lysate

Lane 2: HepG2 whole cell lysate

Lane 3: A549 whole cell lysate


Lane 4: MCF-7 whole cell lysate

Secondary


Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 48 kDa

Observed band size: 48 kDa



Immunohistochemistry of paraffin-embedded human lung cancer using CSB-PA07589A0Rb at dilution of 1:100



Immunohistochemistry of paraffin-embedded human adrenal gland tissue using CSB-PA07589A0Rb at dilution of 1:100