



fla Antibody, HRP conjugated

Product Code	CSB-PA319964LB01BUD
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P11089
Immunogen	Recombinant Borrelia burgdorferi Flagellar filament 41 kDa core protein (1-336AA)
Raised In	Rabbit
Species Reactivity	Borrelia burgdorferi
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Flagellar filament 41 kDa core protein (Flagellin) (41 kDa antigen) (P41), fla
Immunogen Species	Borrelia burgdorferi (strain ATCC 35210 / B31 / CIP 102532 / DSM 4680)
Research Area	Others
Target Names	fla