



INHBA Antibody, HRP conjugated

Product Code	CSB-PA07589B0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P08476
Immunogen	Recombinant Human Inhibin beta A chain protein (311-426AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Inhibin beta A chain (Activin beta-A chain) (Erythroid differentiation protein) (EDF), INHBA
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	INHBA
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.