



CIAPIN1 Antibody

Product Code	CSB-PA738750HA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6FI81
Immunogen	Recombinant Human Anamorsin protein (1-312AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Anamorsin (Cytokine-induced apoptosis inhibitor 1) (Fe-S cluster assembly protein DRE2 homolog), CIAPIN1
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	CIAPIN1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.