



Intermediate capsid protein VP6 Antibody, HRP conjugated

Product Code	CSB-PA541481LB01ROG
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	B1NKS8
Immunogen	Recombinant Rotavirus A Intermediate capsid protein VP6 protein (1-397AA)
Raised In	Rabbit
Species Reactivity	Rotavirus A
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Intermediate capsid protein VP6
Immunogen Species	Rotavirus A
Research Area	Others
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.