




TPO Antibody

Product Code	CSB-PA12279A0RB
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P07202
Immunogen	Recombinant Human Thyroid peroxidase protein (19-846AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Thyroid peroxidase (TPO) (EC 1.11.1.8), TPO
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	TPO

Image



Immunohistochemistry of paraffin-embedded human thyroid tissue using CSB-PA12279A0RB at dilution of 1:100