



Kng1 Antibody

Product Code	CSB-PA012479LA01RA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P08934
Immunogen	Recombinant Rat Kininogen-1 protein (2-110AA)
Raised In	Rabbit
Species Reactivity	Rat
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Kininogen-1 [Cleaved into: Kininogen-1 heavy chain; Bradykinin; Kininogen-1 light chain], Kng1, Kng
Immunogen Species	Rattus norvegicus (Rat)
Research Area	Others
Target Names	Kng1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.