



# Mmp2 Antibody

<b>Product Code</b>	CSB-PA014666LA01RA
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P33436
<b>Immunogen</b>	Recombinant Rat 72 kDa type IV collagenase protein (110-662AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Rat
<b>Tested Applications</b>	ELISA
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	72 kDa type IV collagenase (EC 3.4.24.24) (72 kDa gelatinase) (Gelatinase A) (Matrix metalloproteinase-2) (MMP-2) [Cleaved into: PEX], Mmp2
<b>Immunogen Species</b>	Rattus norvegicus (Rat)
<b>Research Area</b>	Cardiovascular
<b>Target Names</b>	Mmp2
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.