




S100a8 Antibody

Product Code	CSB-PA020641LA01RA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P50115
Immunogen	Recombinant Rat Protein S100-A8 protein (2-89AA)
Raised In	Rabbit
Species Reactivity	Rat
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein S100-A8 (Calgranulin-A) (Migration inhibitory factor-related protein 8) (MRP-8) (p8) (S100 calcium-binding protein A8), S100a8, Mrp8
Immunogen Species	Rattus norvegicus (Rat)
Research Area	Immunology
Target Names	S100a8

Image



Western Blot

Positive WB detected in: Rat spleen tissue, Rat lung tissue

All lanes: S100a8 antibody at 3.2µg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 11 kDa

Observed band size: 11 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.