



Klk7 Antibody

Product Code	CSB-PA012458LA01RA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P36373
Immunogen	Recombinant Rat Glandular kallikrein-7, submandibular/renal protein (25-261AA)
Raised In	Rabbit
Species Reactivity	Rat
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Glandular kallikrein-7, submandibular/renal (rGK-7) (EC 3.4.21.35) (Esterase B) (Kallikrein-related protein K1) (Proteinase A) (RSKG-7) (Tissue kallikrein), Klk7, Klk-7
Immunogen Species	Rattus norvegicus (Rat)
Research Area	Signal Transduction
Target Names	Klk7