



FLT3LG Antibody, Biotin conjugated

Product Code	CSB-PA008734LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P49771
Immunogen	Recombinant Human Fms-related tyrosine kinase 3 ligand protein (27-184AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Fms-related tyrosine kinase 3 ligand (Flt3 ligand) (Flt3L) (SL cytokine), FLT3LG
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	FLT3LG