



CHRM3 Antibody, FITC conjugated

Product Code	CSB-PA005383LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P20309
Immunogen	Recombinant Human Muscarinic acetylcholine receptor M3 protein (253-492AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Muscarinic acetylcholine receptor M3, CHRM3
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	CHRM3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.