



APOA4 Antibody, HRP conjugated

Product Code	CSB-PA001916HB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P06727
Immunogen	Recombinant Human Apolipoprotein A-IV protein (21-396AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Apolipoprotein A-IV (Apo-AIV) (ApoA-IV) (Apolipoprotein A4), APOA4
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	APOA4