



RTN4R Antibody

Product Code	CSB-PA880152ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BZR6
Immunogen	Recombinant Human Reticulon-4 receptor protein (244-473AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Reticulon-4 receptor (Nogo receptor) (NgR) (Nogo-66 receptor), RTN4R, NOGOR
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	RTN4R
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.