




EFTUD2 Antibody

Product Code	CSB-PA617995ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q15029
Immunogen	Recombinant Human 116 kDa U5 small nuclear ribonucleoprotein component protein (1-205AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	116 kDa U5 small nuclear ribonucleoprotein component (Elongation factor Tu GTP-binding domain-containing protein 2) (SNU114 homolog) (hSNU114) (U5 snRNP-specific protein, 116 kDa) (U5-116 kDa), EFTUD2, KIAA0031 SNRP116
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	EFTUD2

Image



Immunohistochemistry of paraffin-embedded human skin tissue using CSB-PA617995ESR1HU at dilution of 1:100

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.