



lecA Antibody, FITC conjugated

Product Code	CSB-PA313124HC01EZX
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q05097
Immunogen	Recombinant Pseudomonas aeruginosa PA-I galactophilic lectin protein (2-122AA)
Raised In	Rabbit
Species Reactivity	Pseudomonas aeruginosa
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	PA-I galactophilic lectin (PA-IL) (Galactose-binding lectin), lecA, pa1L
Immunogen Species	Pseudomonas aeruginosa (strain ATCC 15692 / DSM 22644 / CIP 104116 / JCM 14847 / LMG 12228 / 1C / PRS 101 / PAO1)
Research Area	Others
Target Names	lecA
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.