



# IecA Antibody, Biotin conjugated

<b>Product Code</b>	CSB-PA313124HD01EZX
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q05097
<b>Immunogen</b>	Recombinant Pseudomonas aeruginosa PA-I galactophilic lectin protein (2-122AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Pseudomonas aeruginosa
<b>Tested Applications</b>	ELISA
<b>Form</b>	Liquid
<b>Conjugate</b>	Biotin
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	PA-I galactophilic lectin (PA-IL) (Galactose-binding lectin), IecA, pa1L
<b>Immunogen Species</b>	Pseudomonas aeruginosa (strain ATCC 15692 / DSM 22644 / CIP 104116 / JCM 14847 / LMG 12228 / 1C / PRS 101 / PAO1)
<b>Research Area</b>	Others
<b>Target Names</b>	IecA
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.