



Histone H2B Antibody

Product Code	CSB-PA355820LA01CSJ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P14001
Immunogen	Recombinant Crocodylus niloticus Histone H2B protein (2-93AA)
Raised In	Rabbit
Species Reactivity	Cairina moschata
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Histone H2B
Immunogen Species	Cairina moschata (Muscovy duck)
Research Area	Others