



Capsid protein VP1 Antibody

Product Code	CSB-PA362102LA01HPM
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P07299
Immunogen	Recombinant Human parvovirus B19 Capsid protein VP1 protein (228-781AA)
Raised In	Rabbit
Species Reactivity	Human parvovirus B19
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Human parvovirus B19
Research Area	Others
Target Names	N/A
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.