



S Antibody

Product Code	CSB-PA333052LA01BJO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P25194
Immunogen	Recombinant Bovine coronavirus Spike glycoprotein protein (326-540AA)
Raised In	Rabbit
Species Reactivity	Bovine coronavirus
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Spike glycoprotein (S glycoprotein) (E2) (Peplomer protein) [Cleaved into: Spike protein S1 (90B); Spike protein S2 (90A)], S
Immunogen Species	Bovine coronavirus
Research Area	Others
Target Names	S
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.