



ORF1 Antibody

| | |
|----------------------------|--|
| Product Code | CSB-PA327230LA01HGE |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P33424 |
| Immunogen | Recombinant Hepatitis E virus genotype 1 Non-structural polyprotein pORF1 protein (60-240AA) |
| Raised In | Rabbit |
| Species Reactivity | Hepatitis E virus genotype 1 |
| Tested Applications | ELISA |
| Form | Liquid |
| Conjugate | Non-conjugated |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Non-structural polyprotein pORF1 [Includes: Methyltransferase (EC 2.1.1.-) (EC 2.7.7.-); Putative papain-like cysteine protease (PLP) (EC 3.4.22.-); NTPase/helicase (EC 3.6.4.-); RNA-directed RNA polymerase (RdRp) (EC 2.7.7.48)] |
| Immunogen Species | Hepatitis E virus genotype 1 (isolate Human/Pakistan/Sar-55) (HEV-1) |
| Research Area | Others |
| Target Names | ORF1 |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |