



# ORF1 Antibody, HRP conjugated

<b>Product Code</b>	CSB-PA327230LB01HGE
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P33424
<b>Immunogen</b>	Recombinant Hepatitis E virus genotype 1 Non-structural polyprotein pORF1 protein (60-240AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Hepatitis E virus genotype 1
<b>Tested Applications</b>	ELISA
<b>Form</b>	Liquid
<b>Conjugate</b>	HRP
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Non-structural polyprotein pORF1 [Includes: Methyltransferase (EC 2.1.1.-) (EC 2.7.7.-); Putative papain-like cysteine protease (PLP) (EC 3.4.22.-); NTPase/helicase (EC 3.6.4.-); RNA-directed RNA polymerase (RdRp) (EC 2.7.7.48)]
<b>Immunogen Species</b>	Hepatitis E virus genotype 1 (isolate Human/Pakistan/Sar-55) (HEV-1)
<b>Research Area</b>	Others
<b>Target Names</b>	ORF1
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.