



ORF1 Antibody, Biotin conjugated

Product Code	CSB-PA327230LD01HGE
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P33424
Immunogen	Recombinant Hepatitis E virus genotype 1 Non-structural polyprotein pORF1 protein (60-240AA)
Raised In	Rabbit
Species Reactivity	Hepatitis E virus genotype 1
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Non-structural polyprotein pORF1 [Includes: Methyltransferase (EC 2.1.1.-) (EC 2.7.7.-); Putative papain-like cysteine protease (PLP) (EC 3.4.22.-); NTPase/helicase (EC 3.6.4.-); RNA-directed RNA polymerase (RdRp) (EC 2.7.7.48)]
Immunogen Species	Hepatitis E virus genotype 1 (isolate Human/Pakistan/Sar-55) (HEV-1)
Research Area	Others
Target Names	ORF1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.