



Vegfc Antibody

Product Code	CSB-PA025835LA01MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P97953
Immunogen	Recombinant Mouse Vascular endothelial growth factor C protein (108-223AA)
Raised In	Rabbit
Species Reactivity	Mouse
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Vascular endothelial growth factor C (VEGF-C) (Flt4 ligand) (Flt4-L) (Vascular endothelial growth factor-related protein) (VRP), Vegfc
Immunogen Species	Mus musculus (Mouse)
Research Area	Cardiovascular
Target Names	Vegfc
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.